



Ver. HB230112

# Hieff Clone™ Universal One Step Cloning Kit

## Product description

This kit is further developed based on the principle of efficient and rapid connection of DNA fragments by topoisomerase. Compared with the traditional T4 ligase, it has the following advantages: 1) fast, and the connection reaction can be completed within only 1-5 min. 2) High efficiency, no self-connection, the positive cloning rate is close to 100%, and there is no need to set up blue and white spot screening; 3) The operation is simple, and it takes only 15-20 min from connecting to the coating plate. Ice bath, heat shock and 1 hour resuscitation are omitted in the operation process. 4) It can connect up to 5 kb of products.

## Components

Components No.	Name	10907ES20 (20T)
10907-A	pESI-T vector (30 ng/μL)	20 μL
10907-B	1 kb control insert (40 ng/μL)	5 μL
10907-C	10× Enhancer	20 μL

## Specifications

Product type	Cloning Kit
Antibiotic-resistant Bacteria	Ampicillin (AmpR)
Methods of cloning	TOPO™-TA
Conditions of carriage	Dry ice

## Storage

The product should be stored at -25°C~-15°C for 1 year.

## Instructions

1. Cloning experiment of control DNA fragment

1.1 Prepare the following DNA solution in sterile micro centrifuge tube. Take 10 μL as an example.

Component	Dosage
10 × Enhancer	1 μL
1 kb control insert (40 ng/μL)	1 μL
pESI-T vector (30 ng/μL)	1 μL
ddH <sub>2</sub> O	7 μL

1.2 Mix the above systems. React at room temperature (20-30 °C) for 5 min.

[Note]: The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally,



the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

1.3 The linked products can be transformed directly or stored at -20 °C.

1.4 Total quantity 10 µL add 100 µL competent cells, gently mix and place at room temperature for 5 minutes.

[Note]: a) 5 µL is also acceptable connecting liquid, add 50 µL in competent cells (the added volume shall not exceed 1/10 of the volume of competent cells).

b) Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

1.5 Add 300-500 µL LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

1.6 Take 200 µL bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100 µL. Lightly flick the suspended bacteria and take all the bacterial liquid (coated plate).

## 2. Cloning experiment of general DNA fragments

The inserted fragment is a product containing A tail, which can be amplified by conventional Taq enzyme (YEASEN, Cat#10101-10106), hot start Taq enzyme (YEASEN, Cat#10110), or long fragment DNA polymerase (YEASEN, Cat#10107ES62). If there is no non-specific band and primer dimer, it can be directly connected for ligation reaction. Otherwise, it is recommended to recycle the glue before use.

[Note]: a) PCR products cannot be phosphorylated.

b) If the amplification template is plasmid, the template plasmid will cause false positive in subsequent experiments. Therefore, it is recommended to recover the PCR product and connect it.

2.1 Prepare the connection system according to the following table. Take 10 µL as an example.

Component	Dosage
10 × Enhancer	1 µL
pESI-T vector(30 ng/µL)	1 µL
Insert clip	0.5-8 µL
ddH <sub>2</sub> O	Up to 10 µL

[Note]: a) the reaction system can be adjusted according to the above proportion according to the specific experimental conditions.

2.2 Refer to the following table for the amount of different fragments inserted:

Insert clip size	Recommended dosage
0.1-1 kb	20-50 ng
1-2 kb	50-100 ng
2-5 kb	100-200 ng

2.3 Mix the above systems. React at room temperature (20-30 °C) for 5 min.

[Note]: The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally, the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

2.4 Total quantity 10 µL add 100 µL competent cells, gently mix and place at room temperature for 5 minutes.

[Note]: a) 5 µL is also acceptable connecting liquid, add 50 µL in competent cells (the added volume shall not exceed 1/10 of the volume of competent cells).



b) Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

2.5 Add 300-500 μL LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

[Note]: Generally, when the commercial competent cells do not exceed 2 kb insertion fragments, enough transformants can be obtained after 10 min recovery. If the competent efficiency is low or the insertion fragments are long and there are few transformants, the recovery time can be increased to 30-60 min to obtain more transformants.

2.6 Take 200 μ L bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100 μL. Lightly flick the suspended bacteria and take all the bacterial liquid coated plate).

2.7 Screening and identification of transformants

2.7.1 Colony / bacterial liquid PCR identification;

2.7.2 Plasmid size identification: select monoclonal and identify according to the plasmid size after extracting the plasmid.

2.7.3 Enzyme digestion identification: select appropriate restriction enzymes for identification according to the cloning experimental design.


2.7.4 Sequencing analysis: optional sequencing primer sequences are as follows:

M13F: TGTA AACGACGGCCAGT

M13R: CAGGAAACAGCTATGACC

[Note]: The positive rate of this product is quite high. Generally, the positive cloning rate is close to 100%. As long as the growing colonies are normal (not contaminated miscellaneous bacteria, and the number of transformants is not too small), they are basically positive clones. Therefore, when the inserted fragment does not exceed 2-3 kb, you can directly select 1-2 bacteria for sequencing without identification.

### pESI-T vector Map





```

CTACCGAAGAAAGGCCACCCGCGTGAAGGTGAGCCAGTGAGTTGATTGTGTA AACGACGGCCAGTGTCTGAGGCTCGCTGCAGTCCTG
GATGGCTTCTTTCCGGGTGGCACTTCCACTCGGTCACCTCACTAACACATTTTGTGCCGGTCACAGACTCCGAGCGACGTCAGGAC
HindIII EcoRV EcoRI EcoRI EcoRV NotI PstI
AAGCTTGATATCGAATTCGCGTGTGCCCTT AGGGCGACACGCGAATTCGATATCGCGGCCGCTGCAG
TTCGAACTATAGCTTAAGCGCACAGCGGGA Insert clip TTCCCGCTGTGCGCTTAAGCTATAGCGCCGGCGGACGTC
M13 Reverse primer
TCAATACTGACGATGTCATAGCTGTTTCCTGTCCATAGCAGAAAGTCAAAGCCTCCGACCGGAGGCTTTTGACTTGATCG
AGTTATGACTGCTACCAGTATCGACAAAGGACAGGTATCGTCTTTAGTTTTTCGGAGGCTGGCTCCGAAAACCTGAACTAGC
    
```

**pESI-T vector sequence**

**ORIGIN**

```

1 ctgtaagtgg tggcctaact acggtcac tagaagaaca gtatttgta tctgcgctct
61 gctgaagcca gttacctcg aaaaagagtt gtagtcttt gatccggcaa acaaccacc
121 gctggtagcg gtggttttt ttgttgcaag cagcagatta cgcgagaaa aaaaggatct
181 caagaagatc ctttgatctt ctaccgaaga aaggcccacc cgtgaagggt agccagtgag
241 ttgattgtgt aaaacgacgg ccagtgtctg aggctcgtc cagtcctgaa gcttgatctc
301 gaattcgcgt gtcgccctta agggcgacac gcgaattcga tatcgggcc gcctgcagtc
361 aactagcag atggtcatag ctggttctc tccatagcag aaagtcaaaa gcctccgacc
421 ggaggctttt gacttgatcg gcacgtaaga ggtccaact tcaccataa tgaataaga
481 tcactaccgg gcgtatctt tgagttatcg agatttcag gagctaagga agctaaaatg
541 agtattcaac atttccgtgt cgccttatt ccctttttg cggcattttg ccttctgtt
601 tttgctcacc cagaacgct ggtgaaagta aaagatgctg aagatcagtt gggtgcacga
661 gtgggttaca tcgaactgga tctcaacagc ggtaagatcc ttgagagttt tcgccccgaa
721 gaacgttttc caatgatgag cacttttaa gttctgctat gtggcgcggt attatcccgt
781 attgacgccg ggcaagagca actcggtcgc gcatacact attctcagaa tgacttggtt
841 gactactcac cagtcacaga aaagcatctt acggatggca tgacagtaag agaattatgc
901 agtgctgcca taacatgag tgataacact gcggccaact tacttctgac aacgatcgga
961 ggaccgaagg agctaaccgc tttttgcac aacatggggg atcatgtaac tcgccttgat
1021 cgttggaac cggagctgaa tgaagccata ccaaacgac agcgtgacac cacgatgcct
1081 gtagcaatgg caacaactg gcgcaacta ttaactggcg aactacttac ttagcttcc
1141 cggcaacaat taatagactg gatggaggcg gataaagttg caggaccact tctgcgctcg
1201 gccttcccg ctggctggtt tattgctgat aaatctggag cgggtgagcg tgggtctcgc
1261 ggtatcattg cagcactggg gccagatggt aagccctccc gtatcgtagt tatctacag
1321 acggggagtc aggcaactat ggatgaacga aatagacaga tcgctgagat aggtgcctca
1381 ctgattaagc attgtaatg agggccaaa tgtaatcacc tggctcacct tcgggtgggc
1441 ctttctcgt tctggcggtt ttccatagg ctccgcccc ctgacgagca tcacaaaaat
1501 cgatgctcaa gtcagagggt gcgaaaccg acaggactat aaagatacca ggcgtttccc
1561 cctggaagct cctcgtgcg ctctcctgtt ccgacctgc cgcttaccgg atacctgtcc
1621 gccttttccc cttcggaag cgtggcgctt tctcatagct cacgctgtag gtatctcagt
1681 tcggtgtagg tcgttcgctc caagctgggc tgtgtgacg aacccccgt tcagcccagc
    
```



```
1741 cgctgcgct tatccgtaa ctatgtctt ggtccaacc cgtaagaca cgacttatc  
1801 cactggcag cagccactg taacaggatt agcagagcga ggtatgtagg cggtgctaca  
1861 gagtt//
```

[Note]: The yellow background is the sequence of polyclonal enzyme digestion sites.

## Notes

1. Please wear the necessary PPE, such lab coat and gloves, to ensure your health and safety!