



Ver. HB230112

# Hieff Clone™ Universal One Step Cloning Kit

## Product description

This kit is further developed based on the principle of efficient and rapid connection of DNA fragments by topoisomerase. Compared with the traditional T4 ligase, it has the following advantages: 1) Fast, and the connection reaction can be completed in only 1-5 min. 2) High efficiency, no self-connection, the positive cloning rate is close to 100%, and there is no need to set up blue and white spot screening; 3) The operation is simple, and it takes only 15-20 min from connecting to the coating plate. Ice bath, heat shock and 1 hour resuscitation are omitted in the operation process. 4) It can connect up to 5 kb of products.

## Components

Components No.	Name	10909ES20 (20T)
10909-A	pESI-Blunt vector (30 ng/μL)	20 μL
10909-B	1 kb control insert (40 ng/μL)	5 μL
10909-C	10× Enhancer	20 μL

## Specifications

Product type	Cloning Kit
Antibiotic-resistant Bacteria	Ampicillin (AmpR)
Methods of cloning	TOPO™-Blunt
Conditions of carriage	Dry ice

## Storage

The product should be stored at -25°C~-15°C for 1 year.

## Instructions

1. Cloning experiment of control DNA fragment

1.1 Prepare the following DNA solution in sterile micro centrifuge tube. Take 10 μL as an example.

Component	Dosage
10 × Enhancer	1 μL
1 kb control insert (40 ng/μL)	1 μL
pESI-Blunt vector (30 ng/μL)	1 μL
ddH <sub>2</sub> O	7 μL

1.2 Mix the above systems. React at room temperature (20-30 °C) for 5 min.

[Note]: The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally,



the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

1.3 The linked products can be transformed directly or stored at -20 °C.

1.4 Total quantity 10 µL add 100 µL competent cells, gently mix and place at room temperature for 5 minutes.

[Note]: a) 5 µL is also acceptable connecting liquid, add 50 µL in competent cells (the added volume shall not exceed 1/10 of the volume of competent cells).

b) Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

1.5 Add 300-500 µL LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

1.6 Take 200 µL bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100 µL. Lightly flick the suspended bacteria and take all the bacterial liquid (coated plate).

## 2. Cloning experiment of general DNA fragments

The inserted fragment is a product containing A tail, which can be amplified by conventional Taq enzyme (YEASEN, Cat#10101-10106), hot start Taq enzyme (YEASEN, Cat#10110), or long fragment DNA polymerase (YEASEN, Cat#10107ES62). If there is no non-specific band and primer dimer, it can be directly connected for ligation reaction. Otherwise, it is recommended to recycle the glue before use.

[Note]: a) PCR products cannot be phosphorylated.

b) If the amplification template is plasmid, the template plasmid will cause false positive in subsequent experiments. Therefore, it is recommended to recover the PCR product and connect it.

2.1 Prepare the connection system according to the following table. Take 10 µL as an example.

Component	Dosage
10 × Enhancer	1 µL
pESI-Blunt vector(30 ng/µL)	1 µL
Insert clip	0.5-8 µL
ddH <sub>2</sub> O	Up to 10 µL

[Note]: a) the reaction system can be adjusted according to the above proportion according to the specific experimental conditions.

2.2 Refer to the following table for the amount of different fragments inserted:

Insert clip size	Recommended dosage
0.1-1 kb	20-50 ng
1-2 kb	50-100 ng
2-5 kb	100-200 ng

2.3 Mix the above systems. React at room temperature (20-30 °C) for 5 min.

[Note]: The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally, the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

2.4 Total quantity 10 µL add 100 µL competent cells, gently mix and place at room temperature for 5 minutes.

[Note]: a) 5 µL is also acceptable connecting liquid, add 50 µL in competent cells (the added volume shall not exceed 1/10 of the volume of competent cells).



b) Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

2.5 Add 300-500  $\mu$ L LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

[Note]: Generally, when the commercial competent cells do not exceed 2 kb insertion fragments, enough transformants can be obtained after 10 min recovery. If the competent efficiency is low or the insertion fragments are long and there are few transformants, the recovery time can be increased to 30-60 min to obtain more transformants.

2.6 Take 200  $\mu$  L bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100  $\mu$ L. Lightly flick the suspended bacteria and take all the bacterial liquid coated plate).

2.7 Screening and identification of transformants

2.7.1 Colony / bacterial liquid PCR identification;

2.7.2 Plasmid size identification: select monoclonal and identify according to the plasmid size after extracting the plasmid.

2.7.3 Enzyme digestion identification: select appropriate restriction enzymes for identification according to the cloning experimental design.


2.7.4 Sequencing analysis: optional sequencing primer sequences are as follows:

M13F: TGTAACGACGGCCAGT

M13R: CAGGAACAGCTATGACC

[Note]: The positive rate of this product is quite high. Generally, the positive cloning rate is close to 100%. As long as the growing colonies are normal (not contaminated miscellaneous bacteria, and the number of transformants is not too small), they are basically positive clones. Therefore, when the inserted fragment does not exceed 2-3 kb, you can directly select 1-2 bacteria for sequencing without identification.

## pESI-T vector Map





CTACCGAAGAAAGGCCACCCGTAAGGTGAGCCAGTGAGTTGATTGTGTAAACGACGGCCAGTGTCTGAGGCTCGCTGCAGTCCTG  
 GATGGCTTCTTTCCGGGTGGGCACTTCCACTCGGTCACTCACTAACACATTTTGTGCCGGTCACAGACTCCGAGCGACGTCAGGAC

HindIII EcoRV EcoRI EcoRI EcoRV NotI PstI

AAGCTTGATATCGAATTGCGGTGTCGCCCTT **Insert clip** AAGGGCGACACGCGAATTCGATATCGCGGCCGCTGCAG  
 TTCCAATATAGCTTAAGCGCACAGCGGGAA TTCCCGCTGTGCGCTTAAGCTATAGCGCCGGCGGACGTC

M13 Reverse primer

TCAATACTGACGATGTCATAGCTGTTTCCTGTCCATAGCAGAAAGTCAAAAGCCTCCGACCGGAGGCTTTTGACTTGATCG  
 AGTTATGACTGCTACCAAGTATCGACAAAGGACAGGTATCGTCTTTCAGTTTTCGGAGGCTGGCTCCGAAAAGTGAAGTACG

## pESI-T vector sequence

### ORIGIN

```

1  ctggaagtgg tggcctaact acggctacac tagaagaaca gtatttgta tctgcgctct
61  gctgaagcca gttacctcg aaaaagagtt gtagctctt gatccggcaa acaaccacc
121  gctggtagcg gtggttttt tgtttgcaag cagcagatta cgcgcagaaa aaaaggatct
181  caagaagatc ctttgatttt ctaccgaaga aaggcccacc cgtgaagggt agccagttag
241  ttgattgtgt aaaacgacgg ccagtgtctg aggctcgctg cagtcctgaa gcttgatatc
301  gaattcgctg gtcgccctta agggcgacac gcgaattcga tatcgcgcc gcctgcagtc
361  aatactgacg atggtcatag ctgtttcctg tccatagcag aaagtcaaaa gcctccgacc
421  ggaggctttt gacttgatcg gcacgtaaga ggttccaact ttaccataa tgaaataaga
481  tcactaccgg gcgtatttt tgagttatcg agatttcag gagctaagga agctaaaatg
541  agtattcaac atttccgtgt cgccttatt ccctttttg cggcattttg ccttctgtt
601  tttgctcacc cagaacgct ggtgaaagta aaagatgctg aagatcagtt gggtgcacga
661  gtgggttaca tcgaactgga tctcaacagc ggtaagatcc ttgagagttt tcgccccgaa
721  gaacgttttc caatgatgag cacttttaaa gttctgctat gtggcgcggt attatcccgt
781  attgacgccg ggcaagagca actcggtcgc cgcatatact attctcagaa tgacttggtt
841  gagtactcac cagtcacaga aaagcatctt acggatggca tgacagtaag agaattatgc
901  agtgctgcca taacatgag tgataacact gcggccaact tacttctgac aacgatcgga
961  ggaccgaagg agctaaccgc tttttgcac aacatggggg atcatgtaac tcgccttgat
1021  cgttggaac cggagctgaa tgaagccata ccaaacgacg agcgtgacac cacgatgcct
1081  gtagcaatgg caacaacgtt gcgcaacta ttaactggcg aactacttac ttagcttcc
1141  cggcaacaat taatagactg gatggaggcg gataaagttg caggaccact tctgcgctcg
1201  gcccttcagg ctggctggtt tattgtctga aaatctggag ccggtgagcg tgggtctcgc
1261  ggtatcattg cagcactggg gccagatggt aagccctccc gtatcgtagt tatctacacg
1321  acggggagtc aggcactat ggatgaacga aatagacaga tcgctgagat aggtgcctca
1381  ctgattaagc attggtaatg agggcccaaa tgtaatcacc tggctcacct tcgggtgggc
1441  ctttctcgtg tgctggcgtt ttccatagg ctccgcccc ctgacgagca tcacaaaaat
1501  cgatgctcaa gtcagagggt gcgaaacccg acaggactat aaagatacca ggcgtttccc
1561  cctggaagct cctctgtgcg ctctcctgtt ccgacctgc cgcttaccgg atacctgtcc
1621  gccttttccc cttcggaag cgtggcgctt tctcatagct cagctgtag gtatctcagt
1681  tcggtgtagg tcgttcgctc caagctgggc tgtgtgcacg aacccccgt tcagcccgac
  
```



```
1741 cgctg'gcct tatccggtaa ctatcgtctt gaggccaacc cggtaagaca cgacttatcg
1801 ccactggcag cagccactgg taacaggatt agcagagcga ggtatgtagg cggtgctaca
1861 gagtt//
```

[Note]: The yellow background is the sequence of polyclonal enzyme digestion sites.

## Notes

1. Please wear the necessary PPE, such lab coat and gloves, to ensure your health and safety!