

Hieff Clone[™] Zero Blunt Cloning Kit

10909ES

INSTRUCTIONS FOR USE

Ver. EN20230918

A large, decorative orange wave graphic that spans the bottom of the page. It consists of several overlapping, rounded shapes in various shades of orange, creating a fluid, organic pattern.

CONTENT

PRODUCT DESCRIPTION.....	1
SPECIFICATIONS.....	1
COMPONENTS.....	1
STORAGE.....	1
INSTRUCTIONS.....	1
PESI-BLUNT VECTOR MAP.....	4
NOTES.....	5

Product description

This kit is further developed based on the principle of efficient and rapid connection of DNA fragments by topoisomerase. Compared with the traditional T4 ligase, it has the following advantages: 1) fast, and the connection reaction can be completed within only 1-5 min. 2) High efficiency, no self-connection, the positive cloning rate is close to 100%, and there is no need to set up blue and white spot screening; 3) The operation is simple, and it takes only 15-20 min from connecting to the coating plate. Ice bath, heat shock and 1 hour resuscitation are omitted in the operation process. 4) It can connect up to 5 kb of products.

For rapid cloning of blunt end PCR products; Rapid sequencing of PCR products after cloning (using M13F/M13R primers).

Specifications

Product type	Cloning Kit
Antibiotic-resistant Bacteria	Ampicillin (AmpR)
Methods of cloning	Blunt
Conditions of carriage	Dry ice

Components

Components No.	Name	10909ES20 (20T)
10909-A	pESI-Blunt vector (30 ng/μL)	20 μL
10909-B	1 kb control insert (40 ng/μL)	5 μL
10909-C	10× Enhancer	20 μL

Storage

The product should be stored at -25°C~-15°C for 1 year.

Instructions

1. Cloning experiment of control DNA fragment

1) Prepare the following DNA solution in sterile micro centrifuge tube. Take 10 μL as an example.

Component	Dosage
10 × Enhancer	1 μL
1 kb control insert (40 ng/μL)	1 μL
pESI-Blunt vector (30 ng/μL)	1 μL
ddH ₂ O	7 μL

2) Mix the above systems. React at room temperature (20-30 °C) for 5 min.

*The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally, the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

3) The linked products can be transformed directly or stored at - 20 °C .

4) Total quantity 10 μL add 100 μL competent cells, gently mix and place at room temperature for 5 minutes.

*5 μL is also acceptable connecting liquid, add 50 μL in competent cells (the added volume shall not exceed 1/10 of the volume of

competent cells).

**Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

5) Add 300-500 μL LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

6) Take 200 μL bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4,000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100 μL . Lightly flick the suspended bacteria and take all the bacterial liquid (coated plate).

2. Cloning experiment of general DNA fragments

The inserted fragment is a product containing A tail, which can be amplified by conventional THigh fidelity DNA polymerase(YEASEN, Cat#10164ES). If there is no non-specific band and primer dimer, it can be directly connected for ligation reaction. Otherwise, it is recommended to recycle the glue before use.

*PCR products cannot be phosphorylated.

**If the amplification template is plasmid, the template plasmid will cause false positive in subsequent experiments. Therefore, it is recommended to recover the PCR product and connect it.

1) Prepare the connection system according to the following table. Take 10 μL as an example.

Component	Dosage
10 \times Enhancer	1 μL
pESI-Blunt vector(30 ng/ μL)	1 μL
Insert clip	0.5-8 μL
ddH ₂ O	Up to 10 μL

*the reaction system can be adjusted according to the above proportion according to the specific experimental conditions.

**Refer to the following table for the amount of different fragments inserted:

Insert clip size	Recommended dosage
0.1-1 kb	20-50 ng
1-2 kb	50-100 ng
2-5 kb	100-200 ng

2) Mix the above systems. React at room temperature (20-30 °C) for 5 min.

*The connection reaction cannot be carried out on ice. The reaction time should not exceed 5 minutes. Generally, the ligation reaction can be completed in 1-2 minutes to obtain enough recombinants.

3) Total quantity 10 μL add 100 μL competent cells, gently mix and place at room temperature for 5 minutes.

*5 μL is also acceptable connecting liquid, add 50 μL in competent cells (the added volume shall not exceed 1/10 of the volume of competent cells).

**Generally, commercial competent cells do not need ice bath and heat shock, and enough transformants can be obtained after being placed at room temperature for 5 minutes. If the efficiency of competent cells is low, it can be carried out according to the standard procedure of ice bath heat shock.

4) Add 300-500 μL LB or SOC medium (without antibiotics), shake at 37 °C 180 rpm for 10 min.

*Generally, when the commercial competent cells do not exceed 2 kb insertion fragments, enough transformants can be obtained after 10

min recovery. If the competent efficiency is low or the insertion fragments are long and there are few transformants, the recovery time can be increased to 30-60 min to obtain more transformants.

5) Take 200 μ L bacterial solution coated plate (LB or SOC solid medium containing ampicillin resistance), cultured overnight (if there are few transformants expected, centrifuged at 4,000 rpm for 1 min to obtain more clones, absorbed and discarded part of the supernatant and retained 100 μ L. Lightly flick the suspended bacteria and take all the bacterial liquid coated plate).

6) Screening and identification of transformants

a. Colony / bacterial liquid PCR identification, a single colony is selected into a colony PCR Mix using a sterile gun tip or toothpick (e.g. 2 \times HieffTM Ultra-Rapid HotStart PCR Master Mix (with Dye) Cat#10157ES), adding primers for direct PCR reaction.

b. Plasmid size identification: select monoclonal and identify according to the plasmid size after extracting the plasmid.

c. Enzyme digestion identification: select appropriate restriction enzymes for identification according to the cloning experimental design.


d. Sequencing analysis: optional sequencing primer sequences are as follows:

M13F: TGTAACGACGGCCAGT

M13R: CAGGAAACAGCTATGACC

*The positive rate of this product is quite high. Generally, the positive cloning rate is close to 100%. As long as the growing colonies are normal (not contaminated miscellaneous bacteria, and the number of transformants is not too small), they are basically positive clones. Therefore, when the inserted fragment does not exceed 2-3 kb, you can directly select 1-2 bacteria for sequencing without identification.

pESI-Blunt vector Map



CTACCGAAGAAAGGCCACCCGTGAAGGTGAGCCAGTGAGTTGATTGT GTAAAACGACGGCCAGT GTCTGAGGCTCG CTGCAG TCCTG
 GATGGCTTCTTCCGGGTGGGCACTTCCACTCGGTCACTCAACTAACACATTTTGTGCTGCCGGTACAGACTCCGAGCGACGTCAGGAC
 HindIII EcoRV EcoRI AAGCTT GATATC GAATTC GCGTGTCGCCCTT AAGGGCGACACGCGAATTC GATATC GCGGCCGC CTGCAG
 TTCGAACTATAGCTTAAGCGCACAGCGGAA **插入片段** TTCCCGCTGTGCGCTTAAGCTATAGCGCCGGCGGACGTC
 M13 Forward primer → PstI
 TCAATACTGACGATG GTCATAGCTGTTTCCTG TCCATAGCAGAAAGTCAAAAGCCTCCGACCGGAGGCTTTTGACTTGATCG
 AGTTATGACTGCTACCAGTATCGACAAAGGACAGGTATCGTCTTTCAGTTTTCCGAGGCTGGCTCCGAAAACCTGAACTAGC
 M13 Reverse primer ←

pESI-blunt vector sequence

ORIGIN

```

1 ctgaagtgg tggcctaact acggctacac tagaagaaca gtatttgta tctgcgctct
61 gctgaagcca gttacctcg aaaaagagtt gtagctctt gatccggcaa acaaccacc
121 gctggtagcg gtggTTTT tgttgcaag cagcagatta cgcgcagaaa aaaaggatct
181 caagaagatc cttgatttt ctaccgaaga aaggcccacc cgtgaagggtg agccagtgag
241 ttgattgtg aaaacgacgg ccagtgtctg aggctcgtg cagtctgaa gcttgatctc
301 gaattcgcgt gtcgccctta agggcgacac gcaattcga tatcgcggcc gcctgcagtc
361 aactactgac atggcatag ctgtttcctg tccatagcag aaagtcaaaa gcctccgacc
421 ggaggctttt gacttgatcg gcacgtaaga ggtccaact ttaccataa tgaataaga
481 tcaactaccg gcgtatttt tgagtatcg agatttcag gagctaagga agctaaatg
541 agtattcaac atttccgtgt cgccttatt ccctttttg cggcatttg ccttctgtt
601 tttgctcacc cagaaacgt ggtgaaagta aaagatgctg aagatcagtt gggtgcagca
  
```

661 gtgggttaca tcgaactgga tctcaacagc ggtaagatcc ttgagagttt tcgccccgaa
721 gaacgttttc caatgatgag cacttttaaa gttctgctat gtggcgcggt attatccctg
781 attgacgccg ggcaagagca actcggctgc cgcatacact attctcagaa tgacttggtt
841 gagtactcac cagtcacaga aaagcatctt acggatggca tgacagtaag agaattatgc
901 agtgctgcca taacatgag tgataacact gcggccaact tacttctgac aacgatcgga
961 ggaccgaagg agctaaccgc tttttgac aacatggggg atcatgtaac tcgcttgat
1021 cgttggaac cggagctgaa tgaagccata ccaaacgacg agcgtgacac cacgatgcct
1081 gtagcaatgg caacaacgtt gcgcaaaact ttaactggcg aactacttac tctagcttcc
1141 cggcaacaat taatagactg gatggaggcg gataaagttg caggaccact tctgcgctcg
1201 gcccttcggt ctggctggtt tattgctgat aaatctggag ccggtgagcg tgggtctcgc
1261 ggtatcattg cagcactggg gccagatggt aagccctcc gtatcgtagt tatctacag
1321 acggggagtc aggcaactat ggatgaacga aatagacaga tcgctgagat aggtgcctca
1381 ctgattaagc attgtaatg agggccaaa tgtaatcacc tggctcacct tcgggtgggc
1441 ctttctcgt tgctggcgtt tttccatagg ctccgcccc ctgacgagca tcacaaaaat
1501 cgatgctcaa gtcagaggtg gcgaaaccg acaggactat aaagatacca ggcgtttccc
1561 cctggaagct cctctgctgc ctctctgtt ccgaccctgc cgcttaccgg atacctgtcc
1621 gccttttccc cttcggaag cgtggcgctt tctcatagct cacgctgtag gtatctcagt
1681 tcggtgtagg tcgttcgctc caagctgggc tgtgtgcacg aacccccgt tcagcccagc
1741 cgctgcgct tatccgtaa ctatctctt gactccaacc cgtaagaca cgacttatcg
1801 cactggcag cagccactgg taacaggatt agcagagcga ggtatgtagg cgtgtctaca
1861 gagtt

//

*The yellow background is the sequence of polyclonal enzyme digestion sites.

Notes

1. Materials to be prepared:

1) Self-prepared reagents (only some are listed) :

a. High fidelity enzyme: 2 × Hieff Canace™ AdvanceFast PCR Master Mix (With Dye) (Cat#10164ES) or other equivalent.

b. Colony PCR mix: 2 × Hieff™ Ultra-Rapid HotStart PCR Master Mix (with Dye) (Cat#10157ES) or other equivalent product.

c. Nucleic Acid dye: YeaRed™ Nucleic Acid Gel Stain (10,000 × in Water) (Cat#10202ES) or other equivalent products.

2) Self-supplied instrument consumables (only listed part) : PCR instrument, horizontal electrophoresis tank, glue cutting instrument, EP tube, etc.

2. For your safety and health, please wear lab coats and disposable gloves for operation.

3. This product is for research use ONLY!



To enable success of our customers
Together to make a healthier and brighter world