



# Anti-V5 Affinity Gel

## Product description

Anti-V5 Affinity Gel is prepared by conjugating high quality murine IgG2b monoclonal antibody with Sepharose 4B gel. This product has a high V5 label fusion protein loading capacity (at least 0.5 mg protein/mL gel) and less non-specific binding of miscellaneous proteins. It can be used for immunoprecipitation (IP) and purification of fusion proteins with V5 labels.

Recommended Dosage: 5  $\mu$ L gel per 500  $\mu$ L protein crude body fluid.

## Components

Components No.	Name	20588ES03	20588ES08	20588ES25	20588ES60
20588	Anti-V5 Affinity Gel	1 mL	5 mL	25 mL	100 mL

## Specifications

Clone	10B5
Isotype	Mouse IgG1
Purity	Protein A
Antibody Concentration	7.5 g antibody/L Gel
Application	Protein purification、 IP
Protein binding capacity	$\geq 0.5$ mg protein/mL Gel
Storage buffer	10 mM Na <sub>3</sub> PO <sub>4</sub> , 150 mM NaCl, 50% glycerin, pH7.4, contains 0.02% (w/v) sodium azide

## Storage

The unopened products should be stably stored at -25°C~ -15°C for 1 year. Do not freeze in solution without glycerin!

## Instructions

1. Fully re-suspension of Anti-V5 affinity purification gel to form a uniform solution as far as possible. Transfer 10  $\mu$ L mixture (5  $\mu$ L gel) to a new centrifuge tube.
2. Add 500  $\mu$ L TBS, gently reinsert the Anti-V5 Affinity Gel, centrifuge at 5,000 rpm for 30 s, discard the supernatant, and repeat the above steps 3-4 times.

Note: Remove supernatant as much as possible, but do not inhale gel.

3. 200-1000  $\mu$ L cell lysate was added to the above precipitation. Incubate slowly at 4°C for 2 hours.
4. Centrifuge at 5,000 rpm for 30 s and transfer the supernatant to a new centrifuge tube for storage for later detection.
5. Wash the precipitation with 0.5mL TBS for 3-4 times to complete the supernatant.
6. Elution: The specific elution conditions selected depend on the characteristics of the target protein or downstream application.



- 6.1 The V5 fusion protein was eluted with 100  $\mu$ L V5 polypeptides (150  $\mu$ g/mL) (Cat#20575ES), incubated slowly at 4°C for 30 min, centrifuged at 5,000 rpm for 30 s, and the supernatant was collected.
- 6.2 For SDS-PAGE electrophoresis, 10  $\mu$ L 2 $\times$  protein loading buffer was added and boiled for 5 min. Centrifuge at 5,000 rpm for 30 s and transfer supernatant carefully for SDS-PAGE detection. The SDS in the loading buffer will destroy the antibodies, so the affinity purification gel cannot be reused.

## Notes

1. For your safety and health, please wear lab coats and disposable gloves for operation.
2. For research use only.